

Histology Innovation for a NEW Generation

Coverslipping

SHURClean™ Downdraft Fume Extractor



Adequate Ventilation. Economical & Portable. Protecting Every Life Story In Your Lab.

Each coverslipped tissue specimen has a life story behind it. The crucial element that you're faced with when coverslipping that tissue specimen is sufficient ventilation from the various histological solvents and reagents that are heavier than air. The Downdraft Fume Extractor offers the perfect ventilation system for histologist to perform manual coverslipping applications that utilize Xylene based mounting medium. The microprocessor controlled fume hood offers an **unobstructed work surface, user-friendly control panel, small footprint that is ductless, portable and easy to clean.**

Product Specifications:

Dimensions: WxHxD(in)(cm)	(25x19x19.5) (63.5x48.26x49.53) (in) (cm)
Weight: lbs(kg)	35 lbs; 16kg
Voltage:	100-120 VAC, 50/60Hz, 1.2 Amp and 220-240VAC, 50/60Hz, .8 Amp

How To Order:

SC-DD	SHURClean™ Downdraft Ductless Fume Hood, microprocessor controlled, for xylene mounting medium fume capture. One SCDD-FR filter combo included. Specifications: 100-120VAC, 50/60Hz, 1.2 Amp
SD-DD-220	SHURClean™ Downdraft Ductless Fume Hood, microprocessor controlled, for xylene mounting medium fume capture. One SCDD-FR filter combo included. Specifications: 220-240VAC, 50/60Hz, .8 Amp
SC-DD-FR	Replacement filter combo includes one xylene filter and one particle filter; for SC-DD Downdraft Fume Hoods.

Related Products:

Consumables:	SHURMark® Microscope Slides & Coverglass
Reagents:	SHURMount™ Coverslipping Mounting Media


More Information

PH: 844.643.1129

www.general-data.com/hc

General Data Healthcare helps labs enhance productivity and improve workflow with innovative products, solutions and support that focus on delivering better patient safety & care.

©2014 General Data Healthcare, Inc. All rights reserved. Specifications are current at the time of publication; however, they are subject to change without notification. All trademarks and registered trademarks are property of General Data Healthcare, Inc., Cincinnati, Ohio, USA.



General Data Healthcare, Inc.
Headquarters: 4354 Ferguson Drive, Cincinnati, OH 45245
TBS Facility: 3014 Croasdaile Drive, Durham, NC 27705

Sales Inquiries:
PH: 844.643.1129
Email: talktous@general-data.com
Web: www.general-data.com/hc

Service:
PH: 800.245.5029
Email: customersupport@general-data.com