

# Forceps Warmer™

Portable Heating Unit



- **Three Wells Provide Space for Multiple Forceps**  
No need to wait on individual forceps heating.
- **Heating Eliminates Paraffin Wax Build-Up**  
Wells are heated to 80°C to melt paraffin wax. Drainage allows for runoff and creates easier tissue orientation at embedding.
- **Protective 2-Piece Casing Prevents Wax Contamination**  
This case design is durable, safe, and resistant to paraffin contamination.
- **Angled Wells for Easy Forceps Retrieval**
- **Portable Addition to Any Embedding Station**

**Protect every life story™** in your lab.

# Product Specifications:

Dimensions: WxDxH (in)(cm)	(5x3.25x3) (12.7x 8.25x7.62) (in)(cm)
Weight : lbs(kg)	3(2)
Voltage:	100-120 VAC, 50/60Hz, 0.5A and 220-240VAC, 50/60Hz, 0.25A

Catalog #	Description
FW-120	Forceps Warmer, w/ 3 heated wells; 100-120 VAC, 50/60Hz, 0.5A
FW-220	Forceps Warmer, w/ 3 heated wells; 220-240 VAC, 50/60Hz, 0.25A

## More Information

**PH: 844.643.1129**

**[www.general-data-healthcare.com](http://www.general-data-healthcare.com)**

*General Data Healthcare helps labs enhance productivity and improve workflow with innovative products, solutions and support that focus on delivering better patient safety & care.*

©2018 General Data Healthcare, Inc. All rights reserved. Specifications are current at the time of publication; however, they are subject to change without notification. All trademarks and registered trademarks are property of General Data Healthcare, Inc., Cincinnati, Ohio, USA.



General Data Healthcare, Inc.  
4043 McMann Road, Cincinnati, OH 45245  
844.643.1129

Email: [talktous@general-data.com](mailto:talktous@general-data.com)  
Web: [www.general-data.com/healthcare](http://www.general-data.com/healthcare)

Service:  
PH: 800.245.5029  
Email: [customersupport@general-data.com](mailto:customersupport@general-data.com)