


SHURDry™ SD-II Slide Dryer



(thermometer not included)

■ **Dries Slides Quickly**

The Slide Dryer II is able to dry slide stainer racks within 20 minutes to 1 hour. Operating temperatures range from ambient room temperature to 75 degrees C.

■ **Accommodates Six Slide Racks**

Large drying space allows up to six slide racks to be moved directly from a stainer.

■ **Small Countertop Footprint**

This compact unit also fits next to the cutting station for use before going into the stainer.

■ **Temperature Control System Monitors Heat Levels**

An intuitive temperature control system uses light indicators to monitor levels and avoid overheating.

■ **Easy Cleaning**

Simply open the lid and lift out the grill for quick and easy cleaning.

Protect every life story™ in your lab.

Product Specifications:

Dimensions: WxDxH (in) [cm]	(11 x 10 x 5) [26 x 28 x 13]
Weight: (lbs) [kg]	(11) [5]
Voltage:	100-120 VAC, 50/60Hz, 6.2A and 220-240VAC, 50/60Hz, 3A

Catalog #	Description
SD-II-120	SHURDry Slide Dryer II, forced air, holds 6 conventional slide stainer racks, overheat protection switch, ambient to 75 degrees C. Specifications: 100-120V, 60Hz, 6.2A.
SD-II-220	SHURDry Slide Dryer II, forced air, holds 6 conventional slide stainer racks, overheat protection switch, ambient to 75 degrees C. Specifications: 220-240V, 50/60Hz, 3A.

More Information

844.643.1129
www.general-data.com/hc

General Data Healthcare helps labs enhance productivity and improve workflow with innovative products, solutions and support that focus on delivering better patient safety & care.

©2023 General Data Healthcare, Inc. All rights reserved. Specifications are current at the time of publication; however, they are subject to change without notification. All trademarks and registered trademarks are property of General Data Healthcare, Inc., Cincinnati, Ohio, USA.



General Data Healthcare, Inc.
4043 McMann Road, Cincinnati, OH 45245
844.643.1129

Email: talktous@general-data.com
Web: www.general-data.com/healthcare

Service:
PH: 800.245.5029
Email: customersupport@general-data.com