

# HISTOorientor™

Specimen De-Wrinkle



- **De-Wrinkles & Removes Bubbles on Paraffin Sections**  
The constant warming temperature at 70°C will quickly and easily heat slides to the necessary temperatures to remove wrinkles and bubbles on your sections.
- **Rapid Attachment of Fresh Frozen Sections**  
The HISTOorientor™ can be kept on your cryostat to ensure solid attachment of frozen sections to slides to prevent loosening of the tissue during staining agitation.
- **Requires Limited Laboratory Space**  
Compact size will not take up valuable laboratory space. The 2"x 2" (5.1x5.1cm) pad allows two slides to be placed on the plate at the same time.
- **Light Indicator indicates ON/OFF for Safety**

**Protect every life story™** in your lab.

# Product Specifications:

|                            |   |
|----------------------------|---|
| Dimensions: WxDxH (in)(cm) | (2x2x3) (6x6x8) (in) (cm)                                   |
| Weight : bs(kg)            | 3(2)  |
| Voltage:                   | 100-120 VAC, 50/60Hz, 0.125A and 220-240VAC, 50/60Hz, 0.63A |

| Catalog # | Description  |
|-----------|--|
| TO-120    | HISTO <b>Orientor</b> ™, with 2"x2" heated surface; 100-120 VAC, 50/60Hz, 0.125A |
| TO-220    | HISTO <b>Orientor</b> ™, with 2"x2" heated surface; 220-240 VAC, 50/60Hz, 0.63A  |

## More Information

**PH: 844.643.1129**

**[www.general-data-healthcare.com](http://www.general-data-healthcare.com)**

*General Data Healthcare helps labs enhance productivity and improve workflow with innovative products, solutions and support that focus on delivering better patient safety & care.*

©2018 General Data Healthcare, Inc. All rights reserved. Specifications are current at the time of publication; however, they are subject to change without notification. All trademarks and registered trademarks are property of General Data Healthcare, Inc., Cincinnati, Ohio, USA.



General Data Healthcare, Inc.  
4043 McMann Road, Cincinnati, OH 45245  
844.643.1129

Email: [talktous@general-data.com](mailto:talktous@general-data.com)  
Web: [www.general-data.com/healthcare](http://www.general-data.com/healthcare)

Service:  
PH: 800.245.5029  
Email: [customersupport@general-data.com](mailto:customersupport@general-data.com)