


Protect every life story™ in your lab.

LaserTrack™

AB12 Cassette Printer



■ Reliable, High Resolution Text & Barcodes

The precision laser used in the AB12 cassette printer creates the crisp edges and high resolution, high-contrast alphanumeric text and 2D lines a barcode requires to produce a 99.5% scan rate throughout the tracking process.

■ Twelve Cassette Magazines for High Volume Labs

12 removable cassette magazines hold up to 70 cassettes each, for a total of 840, allowing high volume labs uninterrupted workflow efficiency.

■ Smooth Operation

Cassettes load easily into magazines, and the simple marking mechanism greatly minimizes jamming, as well as providing smooth, quiet operation.

Product Specifications:

Dimensions: Diameter x H	(21x 37) (53.34 x 93.98) (in) (cm)	Environmental	Normal laboratory temperatures/humidity
Weight: lbs, kg	115 lbs; 52.16 Kg	Connectivity	USB/Ethernet/Serial (to PC)
Voltage	115 V 60 Hz 350 VA draw	Barcode Symbologies	Code 128, Code 39, Datamatrix, PDF417 (truncated)
Resolution	200 lpi	Graphics	General scaled/monochrome bitmap only
Speed	10 seconds cassette to cassette (approx. 6 cassettes per minute)	Media Max Load	12 stacks of 70 cassettes each

Capable of interfacing with:

- Cerner
- Epic
- Sunquest
- MEDITECH
- Other LIS/LIMS systems

G1-205-CL12 LaserTrack™ AB12 Laser Cassette Printer

Includes:

- Power Cord
- USB Cable
- User Manual
- 12 magazines
- Cassette Organizer
- Operating Software

More Information

844.643.1129

www.general-data.com/hc

General Data Healthcare helps labs enhance productivity and improve workflow with innovative products, solutions and support that focus on delivering better patient safety & care.

©2018 General Data Healthcare, Inc. All rights reserved. Specifications are current at the time of publication; however, they are subject to change without notification. All trademarks and registered trademarks are property of General Data Healthcare, Inc., Cincinnati, Ohio, USA.



General Data Healthcare, Inc.
4043 McMann Road, Cincinnati, OH 45245
844.643.1129

Email: talktous@general-data.com
Web: www.general-data.com/healthcare

Service:
PH: 800.245.5029
Email: customersupport@general-data.com