

Protect every life story™ in your lab.

# LaserTrack™ JB-Y4 Cassette Printer



- **Reliable, High Resolution Text & Barcodes**  
The precision laser used in the JB-Y4 cassette printer creates the crisp edges and high resolution, high-contrast alphanumeric text and 2D lines a barcode requires to produce a 99.5% scan rate throughout the tracking process.
- **Integrated Cassette Organizer**  
No double-handling of cassettes; the cassette organizer holds up to 100 cassettes for easy organization and delivery to grossing stations.
- **Four Removable Cassette Magazines**  
4 removable cassette magazines hold up to 70 cassettes each, for a total of 280, and can be easily switched with magazines containing different cassette types or colors.
- **Optimized Cassettes**  
LaserTrack cassettes are specifically designed and optimized for use in LaserTrack cassette printers. This design results in easy cassette loading and effective movement within the instrument.
- **Smooth Operation**  
The optical sensor system and simple marking mechanism work in harmony with the optimized cassettes to greatly minimize jamming, as well as providing smooth, quiet operation.

# Product Specifications:

Dimensions: WxDxH	(24 x 15 x 16) (60.96 x 38.1 x 40.64) (in) (cm)	Environmental	Normal laboratory temperatures/humidity
Weight: lbs, kg	88 lbs; 39.9 Kg	Connectivity	USB/Ethernet/Serial (to PC)
Voltage	110/220 VAC 50/60 Hz	Barcode Symbologies	Code 128, Code 39, Datamatrix, PDF417 (truncated)
Resolution	200 lpi	Graphics	General scaled/monochrome bitmap only
Speed	5 seconds cassette to cassette (approx. 12 cassettes per minute)	Media Max Load	4 stacks of 70 cassettes each

Capable of interfacing with:

- Cerner
- Epic
- Sunquest
- MEDITECH
- Other LIS/LIMS systems

## G1-205-JBY4 LaserTrack™ JB-Y4 Laser Cassette Printer

Includes:

- Power Cord
- USB Cable
- User Manual
- 4 magazines
- Cassette Organizer
- Operating Software

### Cassettes

*Tissue Cassettes*

G1-JTP-100(color)	(WE)	(YW)	(PK)	(AQ)	(BE)	(OE)	(GY)	(GN)	(LC)	(PH)	(TN)
-------------------	------	------	------	------	------	------	------	------	------	------	------

*Biopsy Cassettes*

G1-JBP-100(color)	(WE)	(YW)	(PK)	(AQ)	(BE)	(OE)	(GY)	(GN)	(LC)	(PH)	(TN)
-------------------	------	------	------	------	------	------	------	------	------	------	------

Ordering Example: to order Green Tissue Cassettes, the part number would be G1-JTP-100GN

## More Information

**844.643.1129**

**[www.general-data.com/hc](http://www.general-data.com/hc)**

*General Data Healthcare helps labs enhance productivity and improve workflow with innovative products, solutions and support that focus on delivering better patient safety & care.*

©2020 General Data Healthcare, Inc. All rights reserved. Specifications are current at the time of publication; however, they are subject to change without notification. All trademarks and registered trademarks are property of General Data Healthcare, Inc., Cincinnati, Ohio, USA.



General Data Healthcare, Inc.  
4043 McMann Road, Cincinnati, OH 45245  
844.643.1129

Email: [talktous@general-data.com](mailto:talktous@general-data.com)  
Web: [www.general-data.com/healthcare](http://www.general-data.com/healthcare)

Service:  
PH: 800.245.5029  
Email: [customersupport@general-data.com](mailto:customersupport@general-data.com)