


**Protect every life story**<sup>™</sup> in your lab.



## Optimized Cassettes

LaserTrack cassettes are specifically designed and optimized for use in LaserTrack cassette printers. This design results in effective movement within the instrument that reduces jamming and produces optimal scan rates of nearly 100%.

## 1, 2, or 3-Sided Printing Technology

The LaserTrack PH8 is the **ONLY** cassette printer with the ability to print text and barcodes

on one, two, or three sides of the cassette, allowing more information on the cassette. This unit is capable of meeting ever-changing data capture requirements and minimizing errors at the same time.





## Ancillary Consumables Eliminated

The Phaser cassette is all that is required. Hot foil tape, light bulbs, ink cartridges, printheads, and filters are not necessary.



# **Product Specifications:**

Dimensions: LxDxH (14 x 26 x 19) (345 x 650 x 480) (in) (mm) Weight: lbs, kg 51 lbs; 23 Kg

100 - 240V, 50/60Hz, 3A Voltage: The main power supply voltage

fluctuations should not exceed ± 10%.

No. of Magazines

Characters Alphanumeric, linear barcode or data matrix ASCII codes **Network Connection** USB **Electrical Conformity** 2004/108/EEC **EM Compatibility** 2006/95/EEC Low Tension

2006/42/EEC

Devices

#### Capable of interfacing with:

- Cerner MEDITECH
- Epic · Other LIS/LIMS systems
- Sunquest

#### G1-205-CLD08 LaserTrack™ PH8 Laser Cassette Printer

### Includes:

- Power Cord
- USB Cable
- User Manual
- **Operating Software**

## G1-205-CLD08D LaserTrack™ PH8 Duo Laser Cassette Printer Includes: Duo chute cassette management system moves cassettes to right or left for workflow

- efficiency
- Power Cord
- **USB** Cable
- User Manual



	White	Blue	Green	Gray	Pink	Orange	Yellow	Lilac	Beige	Aqua	Red	Ochre
Cassettes*	$\circ$											
Tissue Cassettes												
G1-DTP-100(color)	(WE)	(BE)	(GN)	(GY)	(PK)	(OE)	(YW)	(LC)	(BG)	(AQ)	(RE)	(OC)
Biopsy Cassettes												
G1-DBP-100(color)	(WE)	(BE)	(GN)	(GY)	(PK)	(OE)	(YW)	(LC)	(BG)	(AQ)	(RE)	(OC)

Ordering Example: to order Green Tissue Cassettes, the part number would be G1-DTP-100GN

## **More Information**

844.643.1129 www.general-data.com/hc

General Data Healthcare helps labs enhance productivity and improve workflow with innovative products, solutions and support that focus on delivering better patient safety & care.

©2020 General Data Healthcare, Inc. All rights reserved. Specifications are current at the time of publication; however, they are subject to change without notification. All trademarks and registered trademarks are property of General Data Healthcare, Inc., Cincinnati, Ohio, USA.



General Data Healthcare, Inc. 844.643.1129

Email: talktous@general-data.com Web: www.general-data.com/healthcare

Service:

PH: 800.245.5029

Email: customersupport@general-data.com