

# **Personal ID**® Thermal Barcode Wristbands



## Quick & Easy to Use

Error-free patient identification is one of the most important steps in ensuring patient safety. Just type the information and print. Add a photo for even more safety. Personal ID barcoded patient wristbands help prevent dangerous medical errors by making patient information immediately available for any patient care activity.

#### Reliable & Durable

Personal ID wristbands are extremely durable and will stand up to hand sanitizers, soapy & foam washes, exposure and wear of a typical hospital stay.

## Scans First Time, Every Time

Patient information is easily readable, photos are clear, and barcodes scan every time with Personal ID wristbands. You'll print high contrast, high resolution linear and two-dimensional barcodes that deliver immediate first-time scans.

#### Ultimate in Patient Comfort

Soft, lightweight, and smooth to the patient's skin, with a unique bracelet design in sizes for patients from infant/pediatric to standard adult sizing. They are latex-free and hypo-allergenic, with no skin-irritating agents or chemical additives.

## **Product Specifications:**

#### **Personal ID AC Patient Wristbands**



These wristbands feature a unique pressure-sensitive adhesive closure. Easy to apply, and can be precision-adjusted for a secure, comfortable fit on a wide range of wrist sizes. The tamper-evident closure tab is self-aligning, so no adhesive is exposed to the patient's skin.

- Tamper-evident adhesive closure
  - Durable and water-resistant
- Highly resistant to soap, water, and hand sanitizers
  - Latex-free and hypoallergenic
- "Bracelet" design for added patient comfort

Part Number	Size (W x L)	Qty. Per Roll	Rolls per Carton	Style	Print Method	Applications
L-WB-AC29330	1.545" x 13.25"	80 (1.5" core)	5	Bracelet	Direct Thermal	Adult standard for healthcare facilities
L-WB-ACC293301	1.545" x 13.25"	80 (1.0" core)	5	Bracelet	Direct Thermal	Adult standard for healthcare facilities
L-WB-AC293303	1.545" x 13.25"	250 (3" core)	4	Bracelet	Direct Thermal	Adult standard for healthcare facilities
L-WB-AC19227	1.25" x 9.125"	80 (1.5" core)	5	Bracelet	Direct Thermal	Infant/pediatric standard for healthcare facilities

### **Personal ID SecureSnap Patient Wristbands**



The unique closure provides lightweight comfort and tamperproof security - the wristband cannot be removed without destroying the closure or wristband itself. These wristbands are available in both standard and outpatient configurations for adults, children and infants.

- Tamper-proof comfort snap closure
- Extremely durable; resistant to abrasion, soap, water, alcohol and hand sanitizers
- Print text and barcodes using any desktop thermal printer
- Latex-free and hypoallergenic
- "Bracelet" design for added patient comfort

Part Number	Size (W x L)	Qty. Per Roll	Rolls per Carton	Style	Print Method	Applications
L-WB-SA29295	1.545" x 12"	100 (1.5" core)	5	Bracelet	Direct Thermal	Adult standard for healthcare facilities
L-WB-SA292952	1.545" x 12"	250 (1.5" core)	2	Bracelet	Direct Thermal	Adult standard for healthcare facilities
L-WB-SS28292	1.125" x 11.5"	100 (1.5" core)	5	Strap	Direct Thermal	Adult standard for healthcare facilities
L-WB-OA29295	1.545" x 12"	100 (1.5" core)	5	Bracelet	Direct Thermal	Adult outpatient for healthcare facilities
L-WB-OA292952	1.545" x 12"	250 (1.5" core)	2	Bracelet	Direct Thermal	Adult outpatient for healthcare facilities
L-WB-OP25260	1.545" x 10.625"	100 (1.5" core)	5	Bracelet	Direct Thermal	Pediatric outpatient for healthcare facilities
L-WB-SI18719	3.5" x 11.25"	100 (1.5" core)	5	Bracelet	Direct Thermal	Infant/pediatric outpatient for healthcare facilities



General Data Healthcare, Inc. 4043 McMann Road, Cincinnati, OH 45245 **844.643.1129**  Email: talktous@general-data.com Web: www.general-data.com/healthcare

Service:

PH: 800.245.5029

Email: customersupport@general-data.com