





Protect every life story™ in your lab.

SHURStain™ 2030 Slide Stainer



-  **Compact Footprint for Versatile Placement**
The SHURStain™ is a small benchtop unit measuring only (w/d/h) 47"x17"x14"; ideal for labs with crowded workspaces and shallow counters.
-  **Innovative Software for High-Level Productivity**
Using the Multi-load Mode, this unit can perform up to 5 H&E or 5 Pap processes simultaneously! It also meets requirements for special and user-specified staining protocols.
-  **Easy to Navigate and Operate**
An informative display offers intuitive feedback during staining; and easily recognizable icons help seamlessly guide the operator through adjustments and setup.
-  **Self-Contained Management and Reporting**
The SHURStain™ comes with PC software that enables you to edit all parameters and programs as well as print reports.

Product Specifications:

General Information		No. of programs in the memory	Stores up to 20 programs, each with up to 50 steps
Power Requirements	12V DC - 3A	No. of reagents in the memory	Maximum 52 (32 programmed and 20 user-configurable)
Classification in accordance with IEC 1010	Protective class 1 Pollution degree 2 Overvoltage category II	Immersion Time	1 s - 59 m 59 s per step
Operating Conditions FOR INDOOR USE ONLY	Operating temperature range: 10 - 40°C Relative air humidity: 10 - 80% non-condensing	Number of Stations	20
Dimensions and Weight		Reagent Stations	Maximum 18
Dimensions: WxDxH (in)(mm)	(47.24 x 17.32 x 14.49) (1200 x 440 x 368) (in)	Reagent Tray Volume	300 ml
Standalone Weight: lbs (kg) Total weight with packaging:	121.25 lbs; 55kg 242.51 lbs; 110kg	Cleaning Stations	Maximum 3
Operating Capacity		Loading Stations	Maximum 2
Processing capacity	Up to 5 racks at a time, depending on the programs, load frequency and device configuration. Simultaneous performance of up to 5 different staining protocols	Unloading Stations	Maximum 3
Load capacity per rack	30 slides	Fume Extraction	Active carbon filter
		Battery Life	2 hours (lithium-ion batteries)

How To Order:

SS-2030 SHURStain™ 2030 Linear Slide Stainer; Specifications: 100/120VAC, 50/60Hz, 3A

Related Products:

Consumables: SHURMark® Microscope Slides and Cover Glass

Reagents: SHURStain™ Routine & Special Stains and Formalin


More Information

844.643.1129

www.general-data.com/hc

General Data Healthcare helps labs enhance productivity and improve workflow with innovative products, solutions and support that focus on delivering better patient safety & care.

©2016 General Data Healthcare, Inc. All rights reserved. Specifications are current at the time of publication; however, they are subject to change without notification. All trademarks and registered trademarks are property of General Data Healthcare, Inc., Cincinnati, Ohio, USA.



General Data Healthcare, Inc.
4043 McMann Road, Cincinnati, OH 45245
844.643.1129

Email: talktous@general-data.com
Web: www.general-data.com/healthcare

Service:
PH: 800.245.5029
Email: customersupport@general-data.com