

# SHURStain™

H&E Stains



■ **Hybrid Hematoxylin for Rapid Staining**

Utilizing properties of both Harris and Gill stains, this hybrid hematoxylin quickly provides crisp, clear nuclear delineation.

■ **Alcoholic-Based Eosin**

The variety of types and selection of SHURStain Eosin provides a wide range of different colors and brightnesses.

■ **No Mucin Staining**

Our SHURStain Clarifying Reagents ensure there will be no masking of cells, so pathologists will have a much easier time providing a diagnosis.

■ **No Metallic Sheen**

Eliminate daily filtering and save time at the staining station.

**Protect every life story™** in your lab.

Catalog #	Description
BB-SHH-1054	SHURStain Hematoxylin, 500 ml, 4/cs
BB-SHH-2857	SHURStain Hematoxylin Plus, 500 ml, 4/cs
BB-SHH-3354	SHURStain Hematoxylin GDH1, 500 ml, 4/cs
BB-SHH-3454	SHURStain Hematoxylin GDH2, 500 ml, 4/cs
BB-SEO-1554	SHURStain Eosin, 500 ml, 4/cs
BB-SAA-1401	SHURStain Alcoholic Clarifying Reagent, 1 Gal, 4/cs
BB-SAC-1201	SHURStain Aqueous Clarifying Reagent, 1 Gal, 4/cs
BB-SSB-1701	SHURStain Bluing Reagent, 1 Gal, 4/cs
BB-SHH-5554	SHURStain Select Hematoxylin, 500 ml, 4/cs.
BB-SSH-5754	SHURStain Select Hematoxylin, EXTRA, 500 ml, 4/cs.
BB-SSE-6054	SHURStain Select Eosin-Y, 500 ml, 4/cs.
BB-SSE-6354	SHURStain Select Eosin LITE, 500 ml, 4/cs.
BB-SSE-6554	SHURStain Select Eosin Trichrome, 500 ml, 4/cs.
BB-SSE-6754	SHURStain Select Eosin-Y Phloxine, 500 ml, 4/cs.
BB-SEA-7154	SHURStain Select Clarity-AQ, Concentrate, 500 ml, 4/cs. or 1 each
BB-SEA-7301	SHURStain Select Clarity-AQ, RTU, 1 Gal, 4/cs
BB-SE7-7654	SHURStain Select Clarity-AL70, Concentrate, 500 ml, 4/cs. or 1 each
BB-SSB-8154	SHURStain Select TruBlu, Concentrate, 500 ml, 4/cs. or 1 each
BB-SSB-8301	SHURStain Select TruBlu, RTU, 1 Gal, 4/cs
BB-SSE-6454	SHURStain Eosin-Y w/Phloxine B (GPA), 500 ml, 4/cs
BB-SGH-14500	SHURStain Tru Hematoxylin Gills I, 1 Gal, 4/cs
BB-SGH-24530	SHURStain Tru Hematoxylin Gills II, 1 Gal, 4/cs
BB-SGH-34560	SHURStain Tru Hematoxylin Gills III, 1 Gal, 4/cs
BB-SHH-A3566	SHURStain Tru Hematoxylin Harris Acidified, 1 Gal, 4/cs
BB-SEO-3101	SHURStain Tru Alcoholic Eosin-Y, 1 Gal, 4/cs
BB-SEO-3501	SHURStain Tru Aqueous Eosin-Y 1 Gal, 4/cs

## More Information

844.643.1129

[www.general-data.com/hc](http://www.general-data.com/hc)

*General Data Healthcare helps labs enhance productivity and improve workflow with innovative products, solutions and support that focus on delivering better patient safety & care.*

©2018 General Data Healthcare, Inc. All rights reserved. Specifications are current at the time of publication; however, they are subject to change without notification. All trademarks and registered trademarks are property of General Data Healthcare, Inc., Cincinnati, Ohio, USA.



General Data Healthcare, Inc.  
4043 McMann Road, Cincinnati, OH 45245  
844.643.1129

Email: [talktous@general-data.com](mailto:talktous@general-data.com)  
Web: [www.general-data.com/healthcare](http://www.general-data.com/healthcare)

Service:  
PH: 800.245.5029  
Email: [customersupport@general-data.com](mailto:customersupport@general-data.com)