

SHURTrim™

Portable Block Trimmer



- **Portable Trimmer Can Be Placed for Convenience**
Removes excess paraffin wax from the cassette for easier insertion into microtome quick-release clamps.
- **Microprocessor-Controlled Temperature**
Use the digital readout to set the exact temperature you need.
- **Grooved Surface Heats Quickly**
SHURTrim™ grooved top surface heats up to 75°C and melts excess paraffin from embedded tissue cassettes. Large work surface allows users to de-wax multiple blocks at the same time.
- **Compact & Ergonomic Design**
Compact size will fit anywhere, and the sloped top provides optimum comfort for both right- and left-handed users.
- **Drainage Cup Collects Excess Paraffin**

Protect every life story™ in your lab.

Product Specifications:

Dimensions: WxDxH (in)(cm)	6 x9x7) (15 x23x18) (in)(cm)
Weight : bs(kg)	5 (3)
Voltage:	100-120 VAC, 50/60Hz, 2A and 220-240VAC, 50/60Hz, 1A

Catalog #	Description
ST-120	SHURTrim™ Paraffin block dewaxer with drain cup holder. Includes: 5/ea drain cups.
ST-220	SHURTrim™ Paraffin block dewaxer with drain cup holder. Includes: 5/ea drain cups.

More Information

PH: 844.643.1129

www.general-data-healthcare.com

General Data Healthcare helps labs enhance productivity and improve workflow with innovative products, solutions and support that focus on delivering better patient safety & care.

©2020 General Data Healthcare, Inc. All rights reserved. Specifications are current at the time of publication; however, they are subject to change without notification. All trademarks and registered trademarks are property of General Data Healthcare, Inc., Cincinnati, Ohio, USA.



General Data Healthcare, Inc.
4043 McMann Road, Cincinnati, OH 45245
844.643.1129

Email: talktous@general-data.com
Web: www.general-data.com/healthcare

Service:
PH: 800.245.5029
Email: customersupport@general-data.com