

Protect every life story™ in your lab.

StainerShield™ XT

Slide Identification Labels

■ Designed Specifically for Lab Use

StainerShield XT will remain on the slide after exposure to extremes of heat and cold, routine and special histology stains, harsh chemicals and reagents such as xylene, toluene, formalin, and various concentrations of alcohol.

■ Barcode and Readable Information Do Not Fade

Because the barcode and readable information retain their clarity after exposure to all of these stains, chemicals, and temperature extremes, the barcode still scans and the specimen is still positively identified.

■ Does Not Require an Overlay

StainerShield XT requires no overlamine flap for image protection, saving valuable time in preparing slides for the staining process.

■ Tested to Library of Congress Standards

This label undergoes such rigorous testing that it passes Library of Congress standards. It is proven to last up to 20 years in archival, so your lab can be confident of meeting and exceeding CAP and CLIA compliance guidelines.



StainerShield XT Thermal Transfer Labels (For High Heat Printers)

Label Name	Cat.#	Description
SSXT 7/8" x 7/8", Cognitive - High Heat Ribbon	KXT-05240-78110	StainerShield XT kit contains 3 rolls of 7/8" x 7/8" labels, 1,800 labels on each (total 5,400 labels), 1 label across the web, wound in on a 1" core. Kit includes one 2.4" x 140m ribbon on a 0.5" core.
SSXT 7/8" x 7/8", Zebra GX430t - High Heat Ribbon	KXT-05200-78110	StainerShield XT kit contains 3 rolls of 7/8" x 7/8" labels, 1,800 labels on each (total 5,400 labels), 1 label across the web, wound out on a 1" core. Kit includes two 2.5" x 74m ribbons on a 0.5 inch core that is 4 inches wide.
SSXT 7/8" x 7/8", GoDEX RT230i - High Heat Ribbon	KXT-05220-78110	StainerShield XT kit contains 3 rolls of 7/8" x 7/8" labels, 1,800 labels on each (total 5,400 labels), 1 label across the web, wound out on a 1" core. Kit includes two 2.2" x 74m ribbons on a 0.5" core.
SSXT BA Label (7/8"x5/8"), Cognitive - High Heat Ribbon	KXT-05240-75110	StainerShield XT kit contains 3 rolls of 7/8" x 5/8" labels, 1,800 labels on each (total 5,400 labels), 1 label across the web, wound in on a 1" core. Kit includes one 2.4" x 140m ribbon on a 0.5" core.
SSXT BA Label (7/8"x5/8"), Zebra GX430t - High Heat Ribbon	KXT-05200-75110	StainerShield XT kit contains 3 rolls of 7/8" x 5/8" labels, 1,800 labels on each (total 5,400 labels), 1 label across the web, wound in on a 1" core. Kit includes two 2.5"x 74m ribbons on a 0.5 inch core that is 4 inches wide.
SSXT V-Plus (No Flap), Cognitive - High Heat Ribbon	KXT-05240-31110	StainerShield XT kit contains 3 rolls of 31/32" x 7/8" labels, 1,800 labels on each (total 5,400 labels), 1 label across the web, wound in on a 1" core. Kit includes one 2.4" x 140m ribbon on a 0.5" core.
SSXT 15/16" x 15/16", Cognitive - High Heat Ribbon	KXT-05240-56110	StainerShield XT kit contains , 3 rolls of 15/16" x 15/16" labels, 1,800 labels on each (total 5,400 labels), 1 label across the web, wound in on a 1" core. Kit includes one 2.4" x 140m ribbon on a 0.5" core.
SSXT 7/8"x7/8", 4-across, Cognitive Cxi 4" & Zebra GX430t - High Heat Ribbon	KXT-05430-78410	StainerShield XT kit contains 1 roll of 7/8" x 7/8" labels with 7,200 labels on the roll, 4 labels across the web, wound out on a 1" core. Kit includes one 4.3" x 50m ribbon on a 0.5" core.

StainerShield Thermal Transfer Labels (For Lower Heat Printers)

Label Name	Cat.#	Description
SSXT 7/8" x 7/8", Cognitive - Lower Heat Ribbon	KXT-A5240-78110	StainerShield kit contains 3 rolls of 7/8" x 7/8" labels, 1,800 labels on each (total 5,400), 1 label across the web, wound in on a 1" core. Kit comes with one 2.4"x140m ribbon on a 0.5" core.
SSXT 7/8" x 7/8", Zebra Desktop-sized Printers - Lower Heat Ribbon	KXT-05246-78110	StainerShield kit contains 3 rolls of 7/8" x 7/8" labels, 1,800 labels on each (total 5,400 labels), 1 label across the web, wound out on a 1" core. Kit includes two 2.4" x 73m ribbons on a 0.5" core that is 4" wide.
SSXT 7/8"x7/8", Industrial-sized Printers - Lower Heat Ribbon	KXT-A1430-78130	StainerShield kit contains 3 rolls of 7/8" x 7/8" labels, 1,800 labels on each (total 5,400 labels), 1 label across the web, wound out on a 3" core. Kit includes three 4.3"x 164' ribbons on a 1" core.
SSXT BA Label (7/8"x5/8"), Zebra Desktop-sized Printers - Lower Heat Ribbon	KXT-A5240-75110	StainerShield kit contains 3 rolls of 7/8" x 5/8" labels, 1,800 labels on each (total 5,400 labels), 1 label across the web, wound in on a 1" core. Kit includes two 2.4" x 73m ribbons on a 0.5" core that is 4" wide.
SSXT V-Plus (No Flap), Zebra GX430t - Lower Heat Ribbon	KXT-A5240-31110	StainerShield kit contains 3 rolls of 31/32" x 7/8" labels, 1,800 labels on each (total 5,400 labels), 1 label across the web, wound in on a 1" core. Kit includes two 2.4" x 73m ribbons on a 0.5" core that is 4 inches wide.
SSXT VE- Lower Heat Ribbon	KXT-A1434-15941	StainerShield kit contains 3 rolls of 7/8" x 15/16" labels, 1,500 labels on each (total 4,500 labels), wound in on a 1" core. Kit includes three 4.3"x 164' ribbons on a 1" core.
SSXT 7/8"x7/8", 4-across, Industrial-sized Printers - Lower Heat Ribbon	KXT-A1434-78430	StainerShield kit contains 1 roll of 7/8" x 7/8" labels with 7,200 labels on the roll, 4 labels across the web, wound out on a 3" core. Kit includes one 4.3"x 164' ribbon on a 1" core.
SSXT 7/8x7/8, 4-across with Leader, Industrial-sized printers - Lower Heat Ribbon	KXT-A1434-7843A	StainerShield kit contains 1 roll of 7/8" x 7/8" labels, 7,800 labels on roll, 4 labels across the web, wound out on a 3" core with a leader. Kit includes a 4.3"x 164' ribbon on a 1" core.
SSXT 7/8"x7/8", 4-across, Various Printers that take 1" label cores and 1" ribbon cores - Lower Heat Ribbon	KXT-A1434-78410	StainerShield kit contains 1 roll of 7/8" x 7/8" labels with 7,200 labels on the roll, 4 labels across the web, wound out on a 1" core. Kit includes a 4.3"x 164' ribbon on a 1" core.

©2018 General Data Healthcare, Inc. All rights reserved. Specifications are current at the time of publication; however, they are subject to change without notification. All trademarks and registered trademarks are property of General Data Healthcare, Inc., Cincinnati, Ohio, USA.



General Data Healthcare, Inc.
4043 McMann Road, Cincinnati, OH 45245
844.643.1129

Email: talktous@general-data.com
Web: www.general-data.com/healthcare

Service:
PH: 800.245.5029
Email: customersupport@general-data.com