

TMD™ Tissue Marking Dye



■ No Mordant or Fixative Required

TMD enables you to apply and dry in one easy step. Our unique formulation includes fixative with the dye in the same bottle. With TMD, there is no need for a separate fixative to be applied after the dye – which is one less step for the PA's.

■ Quick Drying

TMD requires only 30 – 60 seconds to dry. This advanced formula applies directly, dries quickly, and marks permanently.

■ Consistent Quality

With TMD, you will be assured of consistent quality throughout the entire package, from when the bottle is first opened to when the last of the dye is dispensed.

■ Low Viscosity Assures Coverage

The low viscosity allows for a thin, even coating of pigment, ensuring superior adherence and permanent marking of tissue throughout the entire process.



Product Specifications:

TMD-7	Tissue Marking Dye, 7-Color Kit, 2oz., 1 each
TMD-KIT-FT-7x2	Tissue Marking Dye, 7-Color Kit, 2 oz., 4 kits/cs
TMD-BK	Tissue Marking Dye, Black, 8 oz., 4 ea/cs
TMD-BK-FT-2	Tissue Marking Dye, Black, 2 oz., 4 ea/cs
TMD-BL	Tissue Marking Dye, Blue, 8 oz., 4 ea/cs
TMD-BL-FT-2	Tissue Marking Dye, Blue, 2 oz., 4 ea/cs
TMD-G	Tissue Marking Dye, Green, 8 oz., 4 ea/cs
TMD-G-FT-2	Tissue Marking Dye, Green, 2 oz., 4 ea/cs
TMD-O	Tissue Marking Dye, Orange, 8 oz., 4 ea/cs
TMD-O-FT-2	Tissue Marking Dye, Orange, 2 oz., 4 ea/cs
TMD-P	Tissue Marking Dye, Purple, 8 oz., 4 ea/cs
TMD-P-FT-2	Tissue Marking Dye, Purple, 2 oz., 4 ea/cs
TMD-R	Tissue Marking Dye, Red, 8 oz., 4 ea/cs
TMD-R-FT-2	Tissue Marking Dye, Red, 2 oz., 4 ea/cs
TMD-Y	Tissue Marking Dye, Yellow, 8 oz., 4 ea/cs
TMD-Y-FT-2	Tissue Marking Dye, Yellow, 2 oz., 4 ea/cs

More Information

844.643.1129

www.general-data.com/hc

General Data Healthcare helps labs enhance productivity and improve workflow with innovative products, solutions and support that focus on delivering better patient safety & care.

©2020 General Data Healthcare, Inc. All rights reserved. Specifications are current at the time of publication; however, they are subject to change without notification. All trademarks and registered trademarks are property of General Data Healthcare, Inc., Cincinnati, Ohio, USA.



General Data Healthcare, Inc.
4043 McMann Road, Cincinnati, OH 45245
844.643.1129

Email: talktous@general-data.com
Web: www.general-data.com/healthcare

Service:
PH: 800.245.5029
Email: customersupport@general-data.com